

Max3™ Plus Breast Board

The Max 3 Plus Breast Board is made out of carbon fiber and is specifically designed to treat breast and lung patients. The system features large grid panels and an integrated carbon fiber wedge system. The system also accommodates a large variety of accessories to secure patients during treatment.



The Max3 Plus Breast Board features the following benefits:


- Four integrated elevation carbon fiber wedges at 7°, 12°, 17°, and 23°.
- Large grid panels for increased treatment angles.
- Indexing Backstop Cushion to prevent patients from involuntarily sliding.
- Able to accommodate patients up to 350lbs.
- Compatible with Embrace Head Only VersaFrames.
- Compatible with Silverman head supports.
- Accommodates Butterfly Board, Arm Supports and SecureVac™ Bags



Breast Board Accessories



Clocking Plate Arm Support
RT-5035-60 Left Support
RT-5035-70 Right Support



RT-5043
Calibrated Arm Cups



RT-8335
Head Only VersaFrame

Butterfly Boards



RT-7527

Carbon Fiber Extended Butterfly Board



RT-7526

Carbon Fiber Butterfly Board

*All Butterfly Boards include adapter posts to allow for breast board compatibility.
Acrylic Butterfly Boards also available.*

SecureFit™ Bars

SecureFit Bars are also available in or order to securely attach the Max3 to your treatment table.



RT-7010 Carbon Fiber SecureFit™ Bar Varian Exact



RT-5010A SecureFit™ Bar Varian Exact



RT-5018 Elekta™ Bar, RT-5016 Siemens Bar



RT-5012A SecureFit™ Bar II Varian Exact

Indexable sliding bar allows for off-center treatments.