

OPTEK™ OVERLAY SYSTEM



The OPTEK™ system was developed to provide the highest degree of stability for the head, neck and shoulders. The system is ideal for IMRT and conformal treatments in a 360° arc. Dual-configured for adult and pediatric treatments.




R630
OPTEK™ Carbon Fibre Overlay
Shown here with optional shoulder suppression and hand grips.

USE WITH:





S-Type Head & Shoulder Masks (see page 3)



S-Type Head Masks (see page 4)


OPTEK™ System Features:

- A maximum width of 61cm.
- CT-compatible.
- Fits all S-type masks from Klarity, Civco, Qfix and others.
- Offers pediatric head, neck & shoulder thermoplastic.
- Indexable.
- Carbon fiber construction.
- Excellent transmission factors.



Dual-configured for adult and pediatric treatments

USE WITH:




Pediatric S-Type Head & Shoulder Masks (see pages 3 & 5)

INDEXING BARS



RE-6
For Siemens and Elekta



RE-7
For Varian Exact