

**TABELLA DI ASSOCIAZIONE CODICI SORGENTI – PRINCIPALI TOMOGRAFI PET**  
**ASSOCIATION TABLE SOURCES CODES – MAIN PET SCANNERS**

BRAND	MODEL	TYPE	EZIP CODES	ERS CODES
SIEMENS	ECAT BIOGRAPH	SET	1 x EG-0317-1.2M 2 x HEGL-0080-1M	1 X FNT1010 2 X CBK1020
SIEMENS	ECAT BIOGRAPH PICO	FANTOCCIO	1 x EG-0317-2M	FNT1014
SIEMENS	ECAT BIOGRAPH PICO	LINEARE	2 x HEGL-0080-1M	CBK1020
SIEMENS	ECAT BIOGRAPH PICO	SET	1 x EG-0317-2M 2 x HEGL-0080-1M	1 x FNT1014 2 x CBK1020
SIEMENS	ECAT BIOGRAPH TRUE V BIOGRAPH 6 BIOGRAPH Duo	SET	1 x EG-0318-2M 2 x HEGL-0080-1M 1 X GF-068	1 X FNT1012 2 X CBK1020 1 X TST1031
SIEMENS	BIOGRAPH MCT 40	SET	1 x EG-0318-2.2M 2 x HEGL-0080-1.2M	1 X FNT1012 2 X CBK1020 1 X TST1031
GE	DISCOVERY ST	LINEARE	1 x HEGL-0132	1 x CBK2020
GE	DISCOVERY 600 / 610	FANTOCCIO	1 x VQC-068	1 x FNT2090
GE	DISCOVERY 600 / 610	LINEARE	1 x HEGL-0140	1 x CBK2022
GE	DISCOVERY 600 / 610	SET	1 x VQC-068 1 x HEGL-0140	1 x FNT2090 1 x CBK2022
GE	DISCOVERY 690 / 710	FANTOCCIO	1 x VQC-068	1 x FNT2090
GE	DISCOVERY 690 / 710	LINEARE	1 x HEGL-0139	1 x CBK2021
GE	DISCOVERY 690 / 710	SET	1 x VQC-068 1 x HEGL-0139	1 x FNT2090 1 x CBK2021
GE	DISCOVERY IQ	SET	1 x VQC-068 1 x EG-0320	1 x FNT2090 1 x FNT2010
PHILIPS	GEMINI GXL	SET	1 X GF-0227 6 X MMS05-022-10U	1 X SPT3181 6 X SPT3180
PHILIPS	GEMINI TF/LXL	SET	1 x HEGL-0136 6 X MMS05-022-10U	1 X CBK3120 6 X SPT3180
PHILIPS	CPET ALLEGRO	SPOT	1 x GF-0227	1 X SPT3181