

Klarity[®]

Klarity[®] TRAVELER™ **Patient Transfer System** *with advanced HOVER technology*



Klarity Traveler™ Patient Transfer System

Utilizing our advanced hover technology for easy and efficient patient transfer between CT/MR/PET/Linac/Proton/Brachytherapy, etc.




CT



1.5T and 3.0T MRI



Linac



Proton



Brachytherapy

The universal design of the transfer system greatly improves clinics' operation efficiency and their patient throughput by allowing patient setup on the transfer device in various radiotherapy protocols.



for Brachytherapy



for SBRT



for SRS



for R612-AIO


Klarity Traveler™ for Brachytherapy



Adjustable Leg Supports and Easy-to-Lock Applicator Clamping Device.



Hover technology enables easy and safe patient transfer.



The trolley comes with hooks for hanging disposal bags and others.



Applicator Clamping Device can be used with different afterloaders: It can be fixed in any direction and patient remains at same position during transfer and treatment.



Features and Benefits:

- **Easy**



The built-in air pump improves the device integration.



Height memory for quick adjusting to the preset trolley height.
Manual adjustment available too.

- **Safety**



Emergency brake system.



One Footbrake controls four wheels synchronously.

- **The built-in rechargeable battery enables flexible and continuous use of the transfer system.**
- **High quality patient trolley makes the transfer smooth.**

Klarity[®]



Klarity Medical & Equipment Co., Ltd.

No. 14, 3rd Street Shawan, GETDD

Guangzhou 510730, China

Tel: 86-20-8222 6380

Fax: 86-20-8208 8552

Email: sales@klarity-medical.com

Version No.: KL-CY-20210108-A0