



Klarity

All-In-One

Klarity UNO™ All-In-One



R612 Solo Align™ Full Body System




R605 Duo Align™ Full Body System

ALL-IN-ONE

NEW

Klarity UNO™ All-In-One System

An all-in-one baseplate for S-type masks compatible with small-bore machines



- Head & Neck Solution



- Craniospinal Irradiation
Whole Body Solution



- Head & Thorax Solution



- Hip & Pelvis Solution



- Breast & Pelvic Solution



- SRS Solution



- SBRT Solution



- Prone Breast Solution



- Prone Pelvis Solution

Klarity UNO™


The All-In-One Baseplates & Standard Accessories



The Klarity UNO™ All-In-One Baseplate for the Elekta HexaPOD™ System.




Klarity UNO™ All-In-One MR Compatible Baseplate
Code: R616-SGF



Klarity UNO™

FOR HEAD & NECK SOLUTION

- > For Head, Head & Shoulder Positioning
- > Use with S-type Head / H&S Mask, and Head / H&S AccuCushion® or Vacuum Bag



S-type Head & Shoulder Mask
RLG460ST (LiteGreen™)




Dedicated Headrest for AccuCushion®
R507-T (276 X 424 X 56 mm)



Head & Shoulder AccuCushion®
RG560-SH1 (500 X 520 mm)



Time-Saving Cushion Molding Aid
RSX-K612SA (450 X 478 X 60 mm)



Baseplate
R616-SCF



Trunk Mask
RLG902-2436LA (LiteGreen™)

FOR CRANIOSPINAL IRRADIATION WHOLE BODY SOLUTION



THE HEAD & THORAX SOLUTION


with the use of Overhead Arm Positioner



Head & Thorax Mask
RLG322-S42436 (LiteGreen™)



Overhead Arm Positioner
R616-WCF



R616 Leg Positioner Cushion
R616-L1




R616 Foot Holder
R616-L2



THE HIP & PELVIS SOLUTION



Pelvis Mask
RLG902-2436BA (LiteGreen™)



Klarity UNO™

THE BREAST & PELVIC SOLUTION (TWO-ARMS OVERHEAD SOLUTION WITH WING BOARD)

- > Use with Wingspan™ Baseplate, Mask Fixing Frame and Head & Thorax Mask and Vacuum Bag
- > Ideal for neck, chest, and pelvic regions immobilization
- > Ensures treatment accuracy by avoiding metal artifacts during scanning
- > Easy to use



Mask Fixing Frame
R616-NTCF-3



R616 Wingspan™ Baseplate
R616-NTCF



Head & Thorax Mask
RLG322-S42436 (LiteGreen™)



Body Vacuum Bag:
MR7677-25N


THE SRS SOLUTION

The hand grips, leg positioner cushion, and foot holder work together to provide patients greater relaxation and reliability.

- > Frameless, non-invasive
- > Sub-millimeter accuracy guaranteed



Alternative Frames
(A set of 3 types)




SRS Head Support



Headrests for Klarity SRS Systems
R507-P



SRS Consumables Set, Basic
RBSRS-F01



SRS Consumables Set, Premium
RBSRS-F02



Klarity UNO™ the SRS Solution includes:

- The All-In-One Baseplate (R616-SCF)
- The SRS Package (R616-ACF) and SRS Consumables Set

* R616-ACF consists of SRS Head Support, Headrests(3), Alternative Frames(3)

ALL-IN-ONE

Klarity UNO™

THE SBRT SOLUTION

- > Precise and Reliable
- > User-friendly
- > More Comforts



R616 Wingspan™ Baseplate



Vacuum Bag for Body Support



Vacuum Bag for Knee Support



R616 Leg Positioner Cushion



R616 Foot Holder




Knee Bridge

Belly Bridge



Compression Airbag

The use of airbag as a substitute for paddle to enhance patient comfort.



Indexing Bar



Respiratory Belt (Optional)

Klarity UNO™ the SBRT Solution includes:

- The All-In-One Baseplate (R616-SCF)
- The SBRT Package (R616-BT)

* R616-BT does not include the Respiratory Belt which is optional.

THE PRONE BREAST SOLUTION

- > Designed with patient comfort in mind
- > Positioning accuracy and reproducibility guaranteed
- > Flexible for left or right breast treatment


Klarity UNO™ the Prone Breast Solution includes:

- The All-In-One Baseplate (R616-SCF)
- The Prone Breast System (R616-BE5)

* The R615-BE5 system can be used independently or in combination with the R616-SCF basesplate.




Optional Components:



Thorax Mask
RLG902-2436TL (LiteGreen™)



Breast AccuCushion®
RG570-S4



Prone Breast Vacuum Bag
R7652-28NL-0


THE PRONE PELVIS SOLUTION

- > Designed with patient comfort in mind
- > Positioning accuracy and reproducibility guaranteed
- > The use of mask, vacuum bag or AccuCushion greatly increases positioning accuracy and patient comfort

Klarity UNO™ the Prone Pelvis Solution includes:

- The All-In-One Baseplate (R616-SCF)
- The Prone Pelvis System (R616-AE5)

* The R615-AE5 system can be used independently or in combination with the R616-SCF basesplate.



Optional Components:



Pelvis Mask
RLG902-2436BA (LiteGreen™)



Pelvis AccuCushion®
RG570-S3



Prone Pelvis Vacuum Bag
R7653-20NL-0

Klarity UNO™ ALL-IN-ONE SYSTEM

THE CONSUMABLES AND ACCESSORIES



S-type Head Mask
RLG460SHT (LiteGreen™)



S-type Head & Shoulder Mask
RLG460ST (LiteGreen™)



S-type Head & Thorax Mask
RLG322-S42436 (LiteGreen™)



Head & Thorax Mask
RLG322-S424C-L
 (Cutout on Left side)



Head & Thorax Mask
RLG322-S424C-R
 (Cutout on Right side)



Thorax Mask
RLG902-2436TL (LiteGreen™)



Pelvis Mask
RLG902-2436BA (LiteGreen™)



Trunk Mask
RLG902-2436LA (LiteGreen™)



Head & Shoulder Vacuum Bag
R7202-9NLR
 625 X 582 mm, 9 Liters fill



Body Vacuum Bag
R7569-12NLR-0
 945 X 530 mm, 12 Liters fill



Full Body Vacuum Bag
R7231-15NLR-0
 1347 X 520 mm, 15 Liters fill



Prone Breast Vacuum Bag
R7652-28NL-0
 930 X 810 mm, 28 Liters fill



Prone Pelvis Vacuum Bag

R7653-20NL-0

800 X 780 mm, 20 Liters fill



Vacuum Bag for Body support

R7641-25NLB5

1200 X 800 mm , 25 Liters fill



Vacuum Bag for Knee Support

R7136-4NLT-0

800 X 280 mm , 4 Liters fill



Prone Breast Vacuum Bag

MR7681-20NL-0

900 X 850 mm, 20 Liters fill



Body Vacuum Bags

MR7677-25N



Head AccuCushion®

RG550-H2 (405 X 258 mm)



Head AccuCushion®

RG560-SH1 (500 X 520 mm)



Pelvis AccuCushions®

RG570-S3 (500 X 600 mm)



Breast AccuCushions®

RG570-S4 (500 X 600 mm)



Body AccuCushion®

RG570-S (600 X 500 mm)



Dedicated Headrest for AccuCushions®

R507-T (276 X 424 X 56 mm)



Time-Saving Cushion Molding Aid

RSX-K612SA (450 X 478 X 60 mm)



Mask Fixing Frame

R616-NTCF-3



Cushion Indexer

RSX-TZ-616A



Respiratory Belt

R616-AIR