


The  
**XPERT<sup>®</sup> 20**  
Core Specimen Radiography System

**Compact  
Powerful  
Reliable**



**KUBTEC<sup>®</sup>**  
MEDICAL



## The XPERT® 20 Core Specimen Radiography System



Captures the most cores in a single image with minimal bench space, freeing up valuable biopsy space with 20 years of proven reliability.

- **New Highlight Tool** for microcalcifications
- Auto Calibration
- Single click PACS integration (wired or wireless)
- Compatible with any core biopsy sample
- Large imaging area with small footprint
- Optional cart available

## The KUBTEC® Core Biopsy Specimen Imaging Portfolio

KUBTEC Core Specimen Radiography plays a critical role in breast biopsy workflows, significantly enhancing patient experience by prioritizing speed, accuracy, and comfort.



**The XPERT® 20 System**

Small footprint with a large imaging area.



**The XPERT® 40 System**

Specimen X-ray with a toolkit of innovative, best-in-class features designed to help you be more efficient.

<b>Imaging Area</b>	3" x 5" (7 x 12 cm)	5" x 6" (12 x 15 cm)
<b>Max Magnification</b>	1.4x	2x
<b>High Resolution</b>	✓	✓
<b>HD Optical Camera</b>	✗	✓
<b>Benchtop Design</b>	✓	✗
<b>Use Any Container</b>	✓	✓

## Core Specimen Images

